

RDNZ Customer Bulletin

DATE: February 28th 2008

REFERENCE: CB 08/06

Improved Precision of Elecsys® Myoglobin

Current Situation:

As a result of an adjustment of Anti-myoglobin-Ab-biotin concentration in the R1 reagent, an improvement in precision and recovery has been made with the Elecsys® Myoglobin STAT and Elecsys® Myoglobin 18 min assay.

There is now no HAMA interference, intra-assay precision has dropped to $\leq 3\%$, and inter-assay precision is now $\leq 5\%$. All other specifications remain unchanged.

Actions Required:

There is no change to Cat #'s or packaging. There is no specific action you need to take. The improved formulation will be starting from the following lot #'s:

Elecsys Myoglobin STAT (Cat # 11820788122) **from Lot # 181854, exp 01/09**

Elecsys Myoglobin (Cat # 12178214122) **from Lot # 181855, exp 01/09**

These lot #'s are not being released from Germany until March, so may take some time to filter through.

Please contact us on Techline if you would like to discuss this further, or require any additional information. .