

# RDNZ Customer Bulletin

**DATE:** June 6th 2008

**REFERENCE:** CB 08/22

## Cfas Magnesium Urine Set Point for cobas c501

### **Current Situation:**

We have been informed that there is an incorrect conversion factor for calculating the Cfas set point for magnesium in urine samples on the cobas c501 (ACN I695 and 771). This has lead to a 15% under recovery of urinary magnesium. .

### **Actions Required:**

Please download the updated version of the Cfas value sheet if you run the urine magnesium assay.

We apologise for any inconvenience. Please contact us on Techline if you would like to discuss this further, or require any additional information. .

Roche Ref: PB2008/032